

PharmaTainer™ FAQ

PharmaTainer™ Capacities & Headroom

<http://biofluidfocus.com/downloads>

All PharmaTainer bottles and carboys have headroom, a volume of capacity above the nominal rated and calibrated capacity. This space allows for addition of supplementary reagents and mixing.

Volumes for a given size of container are essentially the same, regardless of the material of construction.

Cat. No.	Nominal Volume	Brim Capacity	Headroom %
x00125-05	125 mL	190 mL	52%
x00500-05	500 mL	636 mL	27%
x01000-05	1000 mL	1200 mL	20%
x02000-05	2000 mL	2270 mL	14%
x05000-01	5 L	5.9 L	18%
x10000-01	10 L	11.6 L	16%
x20000-01	20 L	22.8 L	14%



8207 118th Avenue
Largo FL 33773 USA

Tel: 1-727-437-7777
Email: info@biofluidfocus.com

Manufactured by Cellon SA, Bascharage, Luxembourg
™ PharmaTainer is a trademark of Cellon SA