

## Sampling Assembly Specifications

Expected annula volume of assembly units:

### Specifications for container Assemblies

Fill out the below form as much as possible. If you specify an item without a brand, we will suggest one from our approved vendor list. This works well for silicone tubing, hose clamps and vent filters.

PharmaTainer 1			
PharmaTainer 2			
Manifold	The assembly will incorporate a manifold(s) appropriate to the number of sample containers.		
Note: PharmaTainer Sampling Assemblies are based on 48mm Ported caps.			
Number of ports*			
Tubing			2-Port cap can have only 1 dip tube. 3-Port cap can have 2 dip tubes.
Dip Tube Port 1			
Dip Tube Port 2			
Exteral Tubing Port 2 (opposite dip tube 2)			
Exteral Tubing Port 1 (opposite input port)			
Vent Tubing (Port 3)			
Vent Filter			
Note: House brand vent filters are 0.22 $\mu$ and radiation stable.			
Port 1 Terminal Fitting			
Port 2 Terminal Fitting			
Hose clamps			
Packing	Standard: 2 layers of packaging, plus poly carton liner and full labeling on all levels.		
Irradiation	(		
Description of Assembly and comments			

--	--

Your Contact Information

Name	
Position	
Institution	
Address	
Cit, State	
Postal Code	
eMail:	
Telephone	

Retun the completed for to [Sales@Biofluidfocus.com](mailto:Sales@Biofluidfocus.com)  
We will acknowledge receipt and contact you when  
we have drawing and a bill of materils.